


VISION 3

Intelligent design, operation and analysis for field-flow fractionation




WYATT
TECHNOLOGY

VISION™

Look ahead, see more

VISION, the brains behind FFF



VISION

FFF-MALS—field-flow fractionation with online multi-angle light scattering detection—is a powerful tool for characterizing all types of

macromolecules and nanoparticles in solution or suspension. VISION software is the intelligent human interface to an FFF-MALS system built on Wyatt Technology's Eclipse™ FFF controller and DAWN® MALS instrument.

VISION comprises two primary modules:

VISION DESIGN™ carries out physics-based FFF simulations for *in silico* method design and optimization. The physics engine also calculates diffusion coefficients and zeta potential from FFF and EAF4 measurements.



VISION DESIGN

VISION RUN™ is comprehensive software for running FFF methods. It seamlessly coordinates the pump and autosampler front-end with the Eclipse, detectors and ASTRA®, and records FFF, electrical and UV signals for diagnostics and analysis in VISION DESIGN.



VISION RUN




VISION interfaces with **ASTRA**, the industry-leading light scattering analysis software for determination of molar mass, size, conformation, composition, particle concentration and more.

Integrated Workflow

From method development to final results

VISION offers a seamless workflow in a few simple steps, starting with developing and optimizing the flow program *in silico*, then running a sample sequence, followed by data processing, refining the separation method, reporting and consolidating the data files in projects.






VISION DESIGN

FFF Method Development

Let your PC do the work

Have a new sample to analyze? VISION DESIGN eliminates the need to spend inordinate amounts of time running a matrix of experiments to develop an FFF method. Just list the estimated particle sizes, select your FFF channel, membrane and spacer, and specify a cross-flow gradient.


VISION DESIGN immediately displays the predicted fractogram using fundamental FFF theory to calculate retention times. It only takes a few minutes to explore the effect of different cross flow rates and timings to improve resolution and find the best method - all from the comfort of your desk.



Refine the method

The final design can be transferred in a click to VISION RUN to test your method.

Your initial guess of sample and channel parameters may have been inaccurate, but have no concern — results of a physical experiment can be fed back into VISION DESIGN to further refine the method. Usually a single test run suffices to give VISION DESIGN the information it needs to help you fully optimize your method.





VISION RUN


FFF Instrument Control

Enter the cockpit

VISION RUN is a sophisticated control center for the Eclipse FFF system. The dashboard shows the system status in real time with direct access to instrument configurations and sample sequences. Status information from the Eclipse instrument's Health Indicators are displayed; should any warnings come up, you will receive actionable recommendations.

Expand your horizons

VISION RUN uses the ICF instrument control framework for Agilent® modules from the 1260 series. The complete feature sets of autosamplers, isocratic or quaternary pumps, diode array detectors, fluorescence detectors and fraction collectors are supported to extract the maximum from FFF experiments.





VISION RUN

Automated Sample Sequences

Fly on autopilot

Once your sample sequence table is defined in VISION RUN, including FFF operation and MALS analysis methods, it's time to take flight. Sit back and trust VISION to run on autopilot, orchestrating the entire FFF-MALS measurement.

Should the need arise, system health indicators will alert you to intervene and prevent wasted runs. Decided to add more samples to a running sequence? Simply create another line in the table and VISION RUN will incorporate the additional samples.

Sequence Table (micelles.S)

Operator and Project

Operator: Standard Operator

Project: Standard Project

Signal Selection: default

Data Output & Method

EMDF Data Path: L:\C

VISION RUN Method Path: CAP

Store ASTRA Data: ☒ same folder as EMDF ☐ separate folder

ASTRA Data Path:

Sequence Table

Run (15% done) 13 1/1 SC_1_2_0.35c 5 min_0.35-0l 15 min_0.5 Vd.M 01:01 h | 06:45 h

	Enabled	Vial	Volume [μL]	Inj./Vial	Sample Name	VISION RUN Method	ASTRA Method	Comet Duration [min]	Detecto Flow	ASTRA Name
1	<input checked="" type="checkbox"/>	11	15.00	1	17001	SC_1_2_0.35c 5 min_0.35-0l 15 min_0...	FFF	0.00	0.50	264_2020-09-14 15-04 - 17001 SC_1_2_0.35c 5 min_0.35...
2	<input checked="" type="checkbox"/>	13	15.00	1	blank	SC_1_2_0.35c 5 min_0.35-0l 15 min_0...	FFF	0.00	0.50	265_2020-09-14 15-04 - blank SC_1_2_0.35c 5 min_0.35...
3	<input checked="" type="checkbox"/>	11	15.00	1	17001	SC_1_2_0.35c 5 min_0.35-0l 15 min_0...	FFF	0.00	0.50	266_2020-09-14 15-04 - 17001 SC_1_2_0.35c 5 min_0.35...
4	<input checked="" type="checkbox"/>	13	15.00	1	blank	SC_1_2_0.35c 5 min_0.35-0l 15 min_0...	FFF	0.00	0.50	267_2020-09-14 15-04 - blank SC_1_2_0.35c 5 min_0.35...
5	<input checked="" type="checkbox"/>	11	15.00	1	17001	SC_1_2_0.35c 5 min_0.35-0l 15 min_0...	FFF	0.00	0.50	268_2020-09-14 15-04 - 17001 SC_1_2_0.35c 5 min_0.35...
6	<input checked="" type="checkbox"/>	13	15.00	1	blank	SC_1_2_0.35c 5 min_0.35-0l 15 min_0...	FFF	0.00	0.50	269_2020-09-14 15-04 - blank SC_1_2_0.35c 5 min_0.35...
7	<input checked="" type="checkbox"/>	12	15.00	1	17-DP-003	SC_1_2_0.35c 5 min_0.35-0l 15 min_0...	FFF	0.00	0.50	270_2020-09-14 15-04 - 17-DP-003 SC_1_2_0.35c 5 min...
8	<input checked="" type="checkbox"/>	13	15.00	1	blank	SC_1_2_0.35c 5 min_0.35-0l 15 min_0...	FFF	0.00	0.50	271_2020-09-14 15-04 - blank SC_1_2_0.35c 5 min_0.35...
9	<input checked="" type="checkbox"/>	12	15.00	1	17-DP-003	SC_1_2_0.35c 5 min_0.35-0l 15 min_0...	FFF	0.00	0.50	272_2020-09-14 15-04 - 17-DP-003 SC_1_2_0.35c 5 min...

Total solvent needed: 460.12 mL

Insert Append Delete Clear Save Sequence Save Sequence As Abort Sequence OK Cancel

Regulatory Compliance

VISION offers an optional 21 CFR Part 11 compliance package, including IQ/OQ documents and procedures.

- Administrator, researcher, technician and guest access levels
- Full audit trails
- Secure, joint SQL database for VISION and ASTRA data
- Local or remote database connectivity
- Data integrity validation
- Full IQ/OQ procedures and documentation validation






ASTRA

ASTRA

Comprehensive characterization of FFF-MALS data



Building on over 40 years of light scattering research and development, ASTRA includes an unrivaled range of analysis features to characterize macromolecules and nanoparticles.

Learn more about ASTRA's capabilities at www.wyatt.com/astra-benefits




VISION DESIGN

FFF Data Analysis

Particle size distributions

Even without light scattering, VISION DESIGN can read in FFF data collected by VISION RUN to determine the hydrodynamic radius (R_h) and its distribution in the sample based on FFF retention time.


Choose from a calculation based on fundamental FFF theory or calibration with reference standards. The detailed, high-resolution R_h distribution complements FFF-MALS, batch dynamic light scattering (DLS) and nanoparticle tracking analysis (NTA).



VISION compares hydrodynamic radius distributions from FFF theory (red) with geometric radius distributions from MALS (green).

Electrophoretic mobility and zeta potential

VISION calculates zeta potential and electrophoretic mobility from EAF4 measurements made with different applied electric fields. Data treatment is fully automated, including baseline and peak selection plus determination of peak retention time and its shift with the electrical field.





Project and Data Management



Left to Right

Geoffrey K. Wyatt, Chief Executive Officer
Dr. Philip J. Wyatt, Chairman of the Board
Clifford D. Wyatt, President

Wyatt Technology provides absolute macromolecular and nanoparticle characterization solutions by developing the finest instrumentation and services to chemical, petrochemical, pharmaceutical, biotechnological and academic laboratories worldwide. We delight our customers with unparalleled levels of service and support, facilitating their cutting-edge research and development efforts.

VISION is one of many tools in Wyatt's Light Scattering Toolkit used to characterize proteins, polymers and nanoparticles in solution.


Learn more at www.wyatt.com

System data traces

An invaluable function of VISION DESIGN is the overlay of system data traces for experiments within a sequence or project. Deviations in flow rate, pressure, or any other data trace are readily identified to validate results and facilitate troubleshooting.

Smart project administration

VISION organizes all FFF and MALS data as projects, for convenient review and replication, as well as comparison and reporting. Projects can be merged and experiments added or deleted at will, making this a powerful and flexible way to handle large sets of FFF experiments.



Learn more about FFF-MALS: www.wyatt.com/fff-mals

Copyright ©2021, Wyatt Technology Corporation. All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form by any means, electronic, mechanical, photocopying, recording, or otherwise, without the prior written permission of Wyatt Technology Corporation.

One or more of Wyatt Technology Corporation's trademarks or service marks may appear in this publication. For a list of Wyatt Technology Corporation's trademarks and service marks, please see <https://www.wyatt.com/about/trademarks>.

© Wyatt Technology | W2000D

